

**Hear from all directions
and stay focused while
actively participating in
conversations in every
listening environment.**



Starkey logo, Starkey, Evolv and Evolv logo are trademarks of Starkey Laboratories, Inc.

evolv^{AI}
CROS/BiCROS System
Made for single-sided
hearing loss



www.Starkey.com

Global Headquarters
6700 Washington Ave. S.
Eden Prairie, MN 55344
1.800.328.8602

©2021 Starkey Laboratories, Inc. All Rights Reserved. 86582-000 8/21 BROC2942-00-EE-ST



Single-sided hearing loss solutions are effortless with Evolv AI.

Our CROS/BiCROS System paired with our most advanced hearing technology and Starkey Sound™ is making hearing effortless.

Evolv AI product family offers a wireless CROS/BiCROS system in select styles. Loaded with our most innovative technology and clearest sound, these hearing aids seamlessly adapt to your lifestyle. The Evolv AI CROS/BiCROS System enhances audibility and speech understanding, ensuring you can hear and experience the world around you.

CROS/BiCROS solutions. Now, available with Evolv AI.

Our always-on, always automatic, and always effortless hearing technology is designed so you won't miss out on the conversations and sounds that bring you joy.



Starkey's easy-to-use CROS/BiCROS System in rechargeable and non-rechargeable styles delivers our clearest Starkey Sound with 55 million personalized adjustments every hour.



CROS

If you're unable to hear in one ear and have normal hearing in the other ear, you may benefit from a CROS solution. You'll be able to fully experience sound in your weaker-hearing ear by way of your healthy ear.

How it Works

Challenging listening environments now have a solution. Our CROS System comes equipped with a microphone that picks up sounds and voices from the side of your weaker-hearing ear, then transmits them to a hearing aid receiver that's fitted on your healthy ear. This innovative solution ensures you won't miss the important moments and sounds in your life.



BiCROS

If you're unable to hear in one of your ears and your other ear has hearing loss, you may benefit from a BiCROS solution. This will allow you to experience sound in your ear with little to no hearing and your better ear as well.

How it Works

Our BiCROS System sends sound from a microphone, placed in your ear with little to no hearing, to a hearing aid that is worn in your better-hearing ear. This sound can be personalized to your hearing loss, allowing you to hear sounds from both sides.

Your hearing care professional can help determine what solution will best fit your needs and lifestyle.